

Dräger Alcotest 6810 med

The Dräger Alcotest 6810 med is a special version of the proven Dräger Alcotest series. It has been designed for diagnostic applications, where the CE marking for medical products is necessary. In addition to the electrochemical DrägerSensors high quality and proven reliability, the Dräger Alcotest 6810 med offers a variety of features that ease daily practices within the medical field.



ST-75-2006

Application areas

The Dräger Alcotest 6810 med has been adapted to determine breath alcohol concentrations especially in medical applications. It may be used, for example, to diagnose vascular injuries (about others in the urology during transurethral resections of the prostate) after previous application of an ethanol containing rinsing solution. It is also used for testing patients for previous alcohol consumption at the emergency reception area or before methadone substitution.

Features

The speed and long-term stability of the electrochemical DrägerSensor leads to a short waiting time for the first measurement after the instrument has been turned on, a fast determination of the measured value and a small standard deviation. The illuminated instrument display with information such as test number, date, time, measurement result and error messages as well as the coloured LED and acoustical alarms guide the operator easily through the use of the instrument.

Sampling

During the automatic measurement, which is activated by the breath sample, the breath flow is monitored and the minimum volume required is adjusted to ensure deep lung air is taken. The manually activated measurement makes measurements possible even with people who are not able to deliver intentionally a breath sample.

Data storage

The integrated data memory can store 250 measurement results. Via the interface the results can be transferred into a PC after the measurement.

CE marking as Medical Device

The Dräger Alcotest 6810 med has a CE marking according to the EU Directive 93/42/EEC for Medical Devices. According to this legislation, this CE marking is necessary for all applications involving "diagnosis, prevention, monitoring, treatment or alleviation of disease". The Dräger Alcotest 6810 med fulfils these requirements for applications such as during operations, at emergency reception areas or during methadone substitution.



ST-1887-2006

Dräger Alcotest 6810 med

TECHNICAL DATA

Measuring principle	Electrochemical sensor
Measurement range	0 to 2.5 mg/L
Sampling conditions (automatic measurement)	At the end of the breath sample
	Breath flow: > 3 L/min
	Breath volume: > 0.3 L
Measurement accuracy	0 to 0.5 mg/L: ± 0.008 mg/L
(standard deviation with an ethanol standard)	> 0.5 mg/L: ± 1.7 % of measured value
Sensitivity drift	Typically 0.6 % of measured value per month
Cross sensitivity to anesthetic gases	None
Calibration interval	6 months
Functional test	(recommended)
Dimensions (H x W x D)	Approx. 142 mm x 75 mm x 36 mm / 5.6" x 2.9" x 1.4"
Weight	Approx. 195 g
Electrical supply	2 x 1.5 V Type Alkaline (Mignon, LR6, AA), approx. 1500 measurements at room temperature, or 2 x 1.2 V NiMH rechargeable batteries (Mignon, LR6, AA)
Interface	RS 232 interface for data transmission to a PC
Classification	Measuring instrument: Class I with measuring function
as per Directive 93/42/EEC, Annex IX	Mouthpieces: Class I
UMDNS-Code	17-475

ORDER INFORMATION

Dräger Alcotest 6810 med (Measuring instrument, 3 mouthpieces, batteries, hand strap and case)	83 19 080
Mouthpieces (100 pieces) with non-return valve	68 11 055
Set rechargeable battery, for use with rechargeable NiMH batteries (special battery cover, 2 rechargeable NiMH batteries)	83 18 565
Power supply for charging NiMH batteries inside of the Dräger Alcotest 6810 med	83 16 991
Ethanol standard solution 1,21 g/L for functional test (500 mL)	67 28 838
AlcoView, PC program for communication with the Dräger Alcotest 6810 med	83 17 915
PC connection cable	83 19 715
Dräger Mobile Printer	83 19 310
Printer paper (5 rolls) 7 year shelf life	83 19 002
Printer communication cable with Mini-USB for use with Dräger Mobile Printer	83 18 657

HEADQUARTERS

Dräger Safety AG & Co. KGaA
Revalstrasse 1
23560 Lübeck, Germany

www.draeger.com

SUBSIDIARIES**AUSTRALIA**

Draeger Safety Pacific Pty. Ltd.
Axxess Corporate Park
Unit 99, 45 Gilby Road
Mt. Waverley Vic 3149
Tel +61 3 92 65 50 00
Fax +61 3 92 65 50 95

CANADA

Draeger Safety Canada Ltd.
7555 Danbro Crescent
Mississauga, Ontario L5N 6P9
Tel +1 905 821 8988
Fax +1 905 821 2565

P. R. CHINA

Beijing Fortune Draeger Safety
Equipment Co., Ltd.
A22 Yu An Rd, B Area,
Tianzhu Airport Industrial Zone,
Shunyi District, Beijing 101300
Tel +86 10 80 49 80 00
Fax +86 10 80 49 80 05

FRANCE

Dräger Safety France SAS
3c route de la Fédération,
BP 80141
67025 Strasbourg Cedex 5
Tel +33 3 88 40 59 29
Fax +33 3 88 40 76 67

ITALY

Draeger Safety Italia S.p.A
Via Longarone 35
20080 Zibido San Giacomo (MI),
Italy
Tel +39 02 90 59 49 1
Fax +39 02 90 00 36 86

MEXICO

Draeger Safety S.A. de C.V.
Av. Peñuelas No. 5 Bodega No. 37
Fraccionamiento Industrial
San Pedrito
Querétaro, Qro México
Tel +52 442 246-1113
Fax +52 442 246-1114

NETHERLANDS

Dräger Safety Nederland B.V.
Edisonstraat 53
2700 AH Zoetermeer
Tel +31 79 344 46 66
Fax +31 79 344 47 90

REP. OF SOUTH AFRICA

Dräger South Africa (Pty) Ltd.
P.O.Box 68601
Bryanston 2021
Tel +27 11 465 99 59
Fax +27 11 465 69 53

SINGAPORE

Draeger Safety Asia Pte Ltd
67 Ayer Rajah Crescent
#06-03
Singapore 139950
Tel +65 68 72 92 88
Fax +65 65 12 19 08

SPAIN

Dräger Safety Hispania S.A.
Calle Xaudaró 5
28034 Madrid
Tel +34 91 728 34 00
Fax +34 91 729 48 99

UNITED KINGDOM

Draeger Safety UK Ltd.
Blyth Riverside Business Park
Blyth, Northumberland
NE24 4RG
Tel +44 1670 352-891
Fax +44 1670 356-266

USA

Draeger Safety, Inc.
101 Technology Drive
Pittsburgh, PA 15275
Tel +1 412 787 8383
Fax +1 412 787 2207